

# Developing lab standard procedures and guidelines to improve uptake of eDNA methods by resource managers: An overview of two national-scale projects (READI-Net and iTrackDNA)

Columbia River Basin AIS Team Meeting  
June 2023

<sup>1</sup>Katy Klymus

<sup>2</sup>Margaret Hunter

<sup>3</sup>Caren C. Helbing

<sup>4</sup>Valerie S. Langlois

<sup>3</sup>Jacob J. Imbery

<sup>4</sup>Julie Couillard

<sup>5</sup>Adam Sepulveda

<sup>1</sup> U.S. Geological Survey, Columbia Environmental Research Center, Columbia, MO, U.S.A.

<sup>2</sup> U.S. Geological Survey, Wetland and Aquatic Research Center, Gainesville, FL, U.S.A.

<sup>3</sup> Department of Biochemistry & Microbiology, University of Victoria, Victoria, British Columbia, Canada

<sup>4</sup> Institut National de la Recherche Scientifique, Quebec City, Quebec, Canada

<sup>5</sup> U.S. Geological Survey, Northern Rocky Mountain Science Center, Bozeman, MT, U.S.A



- **U.S. Geological Survey (USGS)** is a United States federal agency within the Department of the Interior.
- Conducts data collection, monitoring, and scientific understanding of the nation's water, energy, mineral and **biological** resources as well as natural hazards monitoring (i.e. earthquakes, flooding)
- **Provide quality scientific information to the public and policy makers** (<https://www.usgs.gov/>)

# ECALS – eDNA Calibration Study

Invasive Carp Regional Coordinating Committee

**Objective:** Understand eDNA detections of invasive carps in the Chicago Area Waterways System

**Aimed to:**

- Develop new markers
- Understand shedding and degradation rates of carp eDNA
- Assess likelihood of alternative vectors or sources of carp eDNA signal

U.S. Fish and Wildlife Service’s bighead and silver carp eDNA monitoring program

- Monitoring the Great Lakes, Illinois, Upper Mississippi and Ohio river systems
- In conjunction with numerous state agencies
- One of the largest and most well developed eDNA monitoring programs in the U.S.



# USGS labs involved in eDNA research



# Transitioning eDNA science from research to management

- Wildlife and natural resource managers are still hesitant to use eDNA tools:
  - Can eDNA tell me if a species is present in an area?
  - How trustworthy is that data?
- Reasons for limited uptake are many but can be distilled to:
  1. poor communication about the state of the science
  2. lack of trust in the accuracy and reliability of the data
  3. limited integration of decision-makers in the research and development process (e.g., Darling 2015).

# Transitioning eDNA science from research to management

- Wildlife and natural resource managers are still hesitant to use eDNA tools:
  - Can eDNA tell me if a species is present in an area?
  - How trustworthy is that data?
- Reasons for limited uptake are many but can be distilled to:
  1. poor communication about the state of the science
  2. lack of trust in the accuracy and reliability of the data
  3. limited integration of decision-makers in the research and development process (e.g., Darling 2015).
- **Lack of standard methods can lead to uncertainty or concern over quality of results and subsequent interpretation. (#2 above)**



# READI-Net

## Rapid eDNA Assessment and Deployment Initiative & Network

*transitioning eDNA biomonitoring from research to actionable science*

### Technology

Autonomous eDNA samplers



### Implementation

Where & when to sample



Lab analysis standards



Communication



Information framework



# READI-Net

## Rapid eDNA Assessment and Deployment Initiative & Network

*transitioning eDNA biomonitoring from research to actionable science*

### Technology

Autonomous eDNA samplers



### Implementation

Where & when to sample



Lab analysis standards



Communication



Information framework





# READI-Net

## Rapid eDNA Assessment and Deployment Initiative & Network

Goal: Develop standards and guidelines for eDNA analysis workflow among labs that will be participating in the READI-Net network

- Assess repeatability and reproducibility within and across labs via **interlaboratory calibration** studies and the development of within and between lab **proficiency testing**
- Create a **standards and guideline document** for lab analysis of READI-Net samples based off of NAS eDNA database standards

# READI-Net

## Rapid eDNA Assessment and Deployment Initiative & Network



### NAS:

USGS Nonindigenous Aquatic Species (NAS) Database

(<https://nas.er.usgs.gov/>)

Incorporation of eDNA data into this database required the development of standards and guidelines (Ferrante et al. 2022)

- Collection, Sample Processing, Contamination (Collection, Extraction, PCR), PCR Assay (Validation, Optimization, Standard curve), Reporting

# iTrackDNA – Tracking at-risk and invasive species with confidence: Opportunities and challenges of eDNA approaches



**Caren HELBING** Project Shared Leader



**Valerie LANGLOIS** Project Shared Co-Leader




**Jérôme DUPRAS** Project Co-Leader & GE<sup>3</sup>LS Leader




**Louis BERNATCHEZ** Project Co-Leader


iTrackDNA




# iTrackDNA



Tracking at-risk and invasive species with confidence: Opportunities and challenges of eDNA approaches



# iTrackDNA



Tracking at-risk and invasive species with confidence: Opportunities and challenges of eDNA approaches



Support standards development



# iTrackDNA – Interlab calibration studies

**Purpose:** To establish lab protocols and set lab performance baselines

**Phase 1:** Assess variation between labs in running standard curves for qPCR/ dPCR



## Methods:

- Each lab ran multiple replicates of standard qPCR standard curves using the same assay
- All reagents and standards were provided by a central lab
- Variables: labs prepare plates with provided serial dilutions, different qPCR thermocycler
- Assessment: LOD and LOQ of each lab's curve was measured






# iTrackDNA – Interlab calibration studies



iTrackDNA



# iTrackDNA – Interlab calibration studies

- Five-fold serial dilutions ranging from 62,500 to 500 copies/reaction (n=8 technical replicates) and 100 to 0.032 copies/reaction (n=24) were prepared in bulk plus 24 NTCs and aliquoted for distribution to participating labs
- Each lab ran three plates of the serial dilutions and used eLOWQuant (Lesperance et al. 2021) to calculate LOD and LOQ

Plate 1

	1	2	3	4	5	6
A	32,150	6,250	1,250	250	50	50
B	32,150	6,250	1,250	250	50	50
C	32,150	6,250	1,250	250	50	50
D	32,150	6,250	1,250	250	50	50
E	32,150	6,250	1,250	250	50	50
F	32,150	6,250	1,250	250	50	50
G	32,150	6,250	1,250	250	50	50
H	32,150	6,250	1,250	250	50	50


Plate 2

	1	2	3	4	5	6	7	8	9	10	11
A	50	10	10	2	2	2	0.4	0.4	0.4	NTC	NTC
B	50	10	10	2	2	2	0.4	0.4	0.4	NTC	NTC
C	50	10	10	2	2	2	0.4	0.4	0.4	NTC	NTC
D	50	10	10	2	2	2	0.4	0.4	0.4	NTC	NTC
E	50	10	10	2	2	2	0.4	0.4	0.4	NTC	NTC
F	50	10	10	2	2	2	0.4	0.4	0.4	NTC	NTC
G	50	10	10	2	2	2	0.4	0.4	0.4	NTC	NTC
H	50	10	10	2	2	2	0.4	0.4	0.4	NTC	NTC

Plate 3

	1	2	3	4	5	6	7	8
A	0.08	0.08	0.08	0.016	0.016	0.016	NTC	10
B	0.08	0.08	0.08	0.016	0.016	0.016	NTC	10
C	0.08	0.08	0.08	0.016	0.016	0.016	NTC	10
D	0.08	0.08	0.08	0.016	0.016	0.016	NTC	10
E	0.08	0.08	0.08	0.016	0.016	0.016	NTC	10
F	0.08	0.08	0.08	0.016	0.016	0.016	NTC	10
G	0.08	0.08	0.08	0.016	0.016	0.016	NTC	10
H	0.08	0.08	0.08	0.016	0.016	0.016	NTC	10


iTrackDNA




# iTrackDNA – Interlab calibration studies

## Results:

- Total number of participating labs: 17
- USGS READInet program: 7 labs
- Experimental design is consistent with draft Canadian Standards Association (CSA) eDNA standard





First CSA Group  
national eDNA standard



# iTrackDNA – Interlab calibration studies

## LOD and LOQ Results



Statistic	LOD (copies/rxn)	LOD 95% CI Lower	LOD 95% CI Upper	LOQ (copies/rxn)	LOQ 95% CI Lower	LOQ 95% CI Upper
Median±MAD	0.5±0.3	0.4±0.2	0.8±0.6	1.9±1.2	1.3±0.9	2.9±2.0
MAX	2.40	1.90	3.60	9.20	7.00	13.60
MIN	0.10	0.10	0.20	0.40	0.30	0.60

# iTrackDNA – Interlab calibration studies Phase 2

## Methods:


- Multiple assays
- Labs will purchase reagents
- Labs will prepare own standard curves given a specific protocol

**Results:** June 2023

## Future Phases:

- Different enzyme/ mastermixes
- Use of tRNA
- Extraction methods

iTrackDNA



# Conclusions

- Lack of method standardization may contribute to the lack of implementation by managers
- There are a numerous guidelines and standards being produced by various entities (examples: CSA, eDNA Society, USDA Genomics, NAS, Aquanet, Finland, Australia)
- We are utilizing these to develop standards specifically for the READinet and iTRACK programs in order to:
  1. Understand what variables most strongly influence lab variation
    - Phase 1 iTRACK: the 17 labs had very low variation in LOD and LOQ, suggesting good lab technique and instrumentation differences were not problematic
  2. Act as a demonstration to the management community.



# Acknowledgements

## READI-net

Elliot Barnhart  
Sara Eldridge  
Stephen Spear  
Chris Merkes  
Tariq Tajjioui  
Murulee Byappanahalli

Ashley Spoljaric  
Patrick Hutchins  
David Pilliod  
Amanda Boone  
John Stechly  
Caitlin Beaver

Devin Jones  
Nathan Thompson  
Megan Voshage

## iTRACK



# Funding

## READI-net

This project is one among a set of coordinated projects funded by (or in part by) the Bipartisan Infrastructure Law through the US Department of the Interior to advance a nationally coordinated Early Detection and Rapid Response Framework.

## iTRACK



# Questions?

Katy Klymus  
kklymus@usgs.gov